

Inventors

Cheu 09/766,659

May 21, 2004

L3 ANSWER 1 OF 2 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1995:241713 HCAPLUS

DOCUMENT NUMBER: 122:75648

ENTRY DATE: Entered STN: 13 Dec 1994

TITLE: A carbohydrate biosensor surface for the detection of uropathogenic bacteria

AUTHOR(S): Nilsson, Kurt G. I.; Mandenius, Carl-Fredrik

CORPORATE SOURCE: Lund, S-223 70, Swed.

SOURCE: Bio/Technology (1994), 12(13), 1376-8

CODEN: BTCHDA; ISSN: 0733-222X

PUBLISHER: Nature Publishing Co.

DOCUMENT TYPE: Journal

LANGUAGE: English

CLASSIFICATION: 9-1 (Biochemical Methods)

ABSTRACT:

We have developed a new surface for use in biosensors that is based on a gold plate covered with a specific carbohydrate receptor structure. The carbohydrate, Gal $\alpha$ 1-4Gal, was bound covalently via a thioalkylcarboxy-spacer, or adsorbed as a neoglycoprotein, to a two-dimensional gold surface. Both types of surfaces showed high specificity in the binding of the uropathogenic bacteria P-fimbriated Escherichia coli compared to the binding of non-infectious bacteria. The signal to noise ratio is sufficiently high to allow specific detection of the bacteria in biosensor applications.

SUPPL. TERM: Carbohydrate biosensor surface uropathogenic bacteria detection

INDEX TERM: Escherichia coli  
(P-fimbriated; carbohydrate biosensor surface for the detection of uropathogenic bacteria)

INDEX TERM: Albumins, uses  
ROLE: DEV (Device component use); PEP (Physical, engineering or chemical process); PROC (Process); USES (Uses)  
(bovine, carbohydrate biosensor surface for the detection of uropathogenic bacteria)

INDEX TERM: Biosensors  
(carbohydrate biosensor surface for the detection of uropathogenic bacteria)

INDEX TERM: Carbohydrates and Sugars, uses  
ROLE: DEV (Device component use); USES (Uses)  
(carbohydrate biosensor surface for the detection of uropathogenic bacteria)

INDEX TERM: Bacteria  
(uropathogenic, carbohydrate biosensor surface for the detection of uropathogenic bacteria)

INDEX TERM: 7440-57-5, Gold, uses  
ROLE: DEV (Device component use); USES (Uses)  
(carbohydrate biosensor surface for the detection of uropathogenic bacteria)

INDEX TERM: 160294-57-5 160294-57-5D, conjugate with bovine serum albumin  
ROLE: DEV (Device component use); PEP (Physical, engineering or chemical process); PROC (Process); USES (Uses)  
(carbohydrate biosensor surface for the detection of uropathogenic bacteria)

L3 ANSWER 2 OF 2 HCAPLUS COPYRIGHT 2004 ACS on STN  
 ACCESSION NUMBER: 1995:235036 HCAPLUS  
 DOCUMENT NUMBER: 122:4932  
 ENTRY DATE: Entered STN: 10-Dec-1994  
 TITLE: Immobilized carbohydrate biosensor for detection of  
 proteins, viruses, or cells  
 INVENTOR(S): Nilsson, Kurt; Mandenius,  
 Carl-Fredrik  
 PATENT ASSIGNEE(S): Swed.  
 SOURCE: PCT Int. Appl., 14 pp.  
 CODEN: PIXXD2  
 DOCUMENT TYPE: Patent  
 LANGUAGE: English  
 INT. PATENT CLASSIF.:  
 MAIN: G01N033-543  
 SECONDARY: C12Q001-00  
 CLASSIFICATION: 9-1 (Biochemical Methods)  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9424561	A1	19941027	WO 1994-SE343	19940418
W: CA, CZ, JP, RU, US				
RW: AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
EP 648333	A1	19950419	EP 1994-914654	19940418
EP 648333	B1	20020109		
R: DE, ES, FR, GB, IT, SE				
US 6231733	B1	20010515	US 1994-356229	19941219
US 2001017270	A1	20010830	US 2001-766659	20010123
PRIORITY APPLN. INFO.:		SE 1993-1270	A	19930419
		WO 1994-SE343	W	19940418
		US 1994-356229	A1	19941219

# ABSTRACT:

A biosensor is disclosed in which an immobilized carbohydrate or a derivative thereof is used to generate a detectable signal when a protein, virus, or cell is bound to the carbohydrate surface. The sensor is an optical sensor, a piezoelec. sensor, an electrochem. electrode, or a thermistor. A method of binding carbohydrates to a gold surface is also described.

SUPPL. TERM: immobilized carbohydrate biosensor; protein detection  
 immobilized carbohydrate biosensor; virus detection  
 immobilized carbohydrate biosensor; cell detection  
 immobilized carbohydrate biosensor

INDEX TERM: Urinary tract  
 (bacteria; immobilized carbohydrate biosensor for  
 detection of proteins, viruses, or cells)

INDEX TERM: Receptors  
 ROLE: BPR (Biological process); BSU (Biological study,  
 unclassified); BIOL (Biological study); PROC (Process)  
 (carbohydrate sequence; immobilized carbohydrate  
 biosensor for detection of proteins, viruses, or cells)

INDEX TERM: (Biosensors)  
 Cell  
 Electrodes  
 Escherichia coli  
 Immobilization, biochemical

Virus  
(immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Proteins, analysis  
ROLE: ANT (Analyte); ANST (Analytical study)  
(immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Aglycons  
Carbohydrates and Sugars, biological studies  
Glycopeptides  
Glycoproteins, biological studies  
Oligosaccharides  
Thiols, biological studies  
ROLE: DEV (Device component use); THU (Therapeutic use);  
BIOL (Biological study); USES (Uses)  
(immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Hexosamines  
ROLE: DEV (Device component use); THU (Therapeutic use);  
BIOL (Biological study); USES (Uses)  
(residue; immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Bacteria  
(urinary tract; immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Albumins, biological studies  
ROLE: DEV (Device component use); THU (Therapeutic use);  
BIOL (Biological study); USES (Uses)  
(conjugates, with galabiose derivative; immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Amino acids, biological studies  
ROLE: DEV (Device component use); THU (Therapeutic use);  
BIOL (Biological study); USES (Uses)  
(glycosyl, immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Glycoproteins, specific or class  
ROLE: DEV (Device component use); THU (Therapeutic use);  
BIOL (Biological study); USES (Uses)  
(neo-, immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Biosensors  
(optical, immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Biosensors  
(piezoelec., immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: Biosensors  
(thermistor-based, immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: 75281-88-8D, derivs.  
ROLE: DEV (Device component use); RCT (Reactant); THU (Therapeutic use); BIOL (Biological study); RACT (Reactant or reagent); USES (Uses)  
(immobilized carbohydrate biosensor for detection of proteins, viruses, or cells)

INDEX TERM: 7440-57-5, Gold, biological studies 13117-26-5D,

Galabiose, derivs., albumin conjugates 30232-12-3,  
Mercaptopropionic acid  
ROLE: DEV (Device component use); THU (Therapeutic use);  
BIOL (Biological study); USES (Uses)  
(immobilized carbohydrate biosensor for detection of  
proteins, viruses, or cells)


INDEX TERM: 151-51-9, Carbodiimide  
ROLE: RCT (Reactant); RACT (Reactant or reagent)  
(immobilized carbohydrate biosensor for detection of  
proteins, viruses, or cells)

INDEX TERM: 13117-26-5  
ROLE: BPR (Biological process); BSU (Biological study,  
unclassified); BIOL (Biological study); PROC (Process)  
(receptor; immobilized carbohydrate biosensor for  
detection of proteins, viruses, or cells)

INDEX TERM: 50-99-7, D-Glucose, biological studies 50-99-7D,  
D-Glucose, analogs 58-86-6, Xylose, biological studies  
58-86-6D, Xylose, analogs 59-23-4, Galactose, biological  
studies 59-23-4D, Galactose, analogs 131-48-6,  
N-Acetylneuraminic acid 131-48-6D, N-Acetylneuraminic  
acid, analogs 2438-80-4, Fucose 2438-80-4D, Fucose,  
analog 3458-28-4, Mannose 3458-28-4D, Mannose, analogs  
ROLE: DEV (Device component use); THU (Therapeutic use);  
BIOL (Biological study); USES (Uses)  
(residue; immobilized carbohydrate biosensor for  
detection of proteins, viruses, or cells)

L4 ANSWER 1 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN  
 RN 160294-57-5 REGISTRY  
 CN Propanoic acid, 3-[[2-[(4-O- $\alpha$ -D-galactopyranosyl- $\beta$ -D-galactopyranosyl)oxy]ethyl]thio]-, hydrazide (9CI) (CA INDEX NAME)  
 FS STEREOSEARCH  
 MF C17 H32 N2 O12 S  
 SR CA  
 LC STN Files: CA, CAPLUS

Absolute stereochemistry.




\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\*

- 1 REFERENCES IN FILE CA (1907 TO DATE)
- 1 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA
- 1 REFERENCES IN FILE CAPLUS (1907 TO DATE)

L4 ANSWER 2 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN  
 RN 75281-88-8 REGISTRY  
 CN  $\beta$ -D-Galactopyranose, 4-O- $\alpha$ -D-galactopyranosyl- (9CI) (CA INDEX NAME)  
 FS STEREOSEARCH  
 MF C12 H22 O11  
 LC STN Files: BEILSTEIN\*, CA, CAPLUS, TOXCENTER, USPATFULL  
 (\*File contains numerically searchable property data)

Absolute stereochemistry.



\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\*

24 REFERENCES IN FILE CA (1907 TO DATE)

8 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA

24 REFERENCES IN FILE CAPLUS (1907 TO DATE)

L4 ANSWER 3 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN

RN 30232-12-3 REGISTRY

CN Propanoic acid, mercapto- (9CI) (CA INDEX NAME)

OTHER CA INDEX NAMES:

CN Propionic acid, mercapto- (7CI)

OTHER NAMES:


CN Mercaptopropionic acid

MF C3 H6 O2 S

CI IDS, COM

LC STN Files: AGRICOLA, ANABSTR, BIOBUSINESS, BIOSIS, BIOTECHNO, CA, CAOLD, CAPLUS, CASREACT, CIN, EMBASE, IFICDB, IFIPAT, IFIUDB, RTECS\*, TOXCENTER, USPAT2, USPATFULL

(\*File contains numerically searchable property data)



D1-SH

154 REFERENCES IN FILE CA (1907 TO DATE)

28 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA

154 REFERENCES IN FILE CAPLUS (1907 TO DATE)

2 REFERENCES IN FILE CAOLD (PRIOR TO 1967)

L4 ANSWER 4 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN

RN 13117-26-5 REGISTRY

CN D-Galactose, 4-O- $\alpha$ -D-galactopyranosyl- (9CI) (CA INDEX NAME)

OTHER CA INDEX NAMES:

CN Galactose, 4-O- $\alpha$ -D-galactopyranosyl- (7CI)


CN Galactose, 4-O- $\alpha$ -D-galactopyranosyl-, D- (8CI)

OTHER NAMES:

CN  $\alpha$ -D-Galp-(1-4)-D-Gal

CN 4 $\alpha$ -Galactobiose  
 CN Galabiose  
 FS STEREOSEARCH  
 MF C12 H22 O11  
 LC STN Files: AGRICOLA, BEILSTEIN\*, BIOBUSINESS, BIOSIS, CA, CANCERLIT,  
 CAOLD, CAPLUS, MEDLINE, TOXCENTER, USPAT2, USPATFULL  
 (\*File contains numerically searchable property data)

Absolute stereochemistry.



\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\*

99 REFERENCES IN FILE CA (1907 TO DATE)  
 8 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA  
 99 REFERENCES IN FILE CAPLUS (1907 TO DATE)  
 2 REFERENCES IN FILE CAOLD (PRIOR TO 1967)

L4 ANSWER 5 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN

RN 7440-57-5 REGISTRY

CN Gold (8CI, 9CI) (CA INDEX NAME)

OTHER NAMES:

CN A 4631  
 CN A 4953  
 CN AY 5022  
 CN Britecote  
 CN Burnish Gold  
 CN C.I. 77480  
 CN C.I. Pigment Metal 3  
 CN Colloidal gold  
 CN Furuuchi 8560  
 CN G 1402  
 CN Gold 197  
 CN Gold black  
 CN Gold element  
 CN Gold Flake  
 CN Gold Leaf  
 CN Gold Powder  
 CN Palegold 5550  
 CN Perfect Gold  
 CN PH 870  
 CN SG 10NK  
 CN Shell Gold  
 CN TR 1306  
 DR 33019-35-1  
 MF Au

CI COM  
LC STN Files: ADISNEWS, AGRICOLA, ANABSTR, AQUIRE, BIOBUSINESS, BIOSIS,  
BIOTECHNO, CA, CABA, CANCERLIT, CAOLD, CAPLUS, CASREACT, CBNB, CEN,  
CHEMCATS, CHEMLIST, CIN, CSCHM, CSNB, DDFU, DETHERM\*, DRUGU, EMBASE,  
ENCOMPLIT, ENCOMPLIT2, ENCOMPPAT, ENCOMPPAT2, HSDB\*, IFICDB, IFIPAT,  
IFIUDB, IPA, MEDLINE, MRCK\*, MSDS-OHS, NIOSHTIC, PIRA, PROMT, RTECS\*,  
TOXCENTER, ULIDAT, USPAT2, USPATFULL, VTB  
(\*File contains numerically searchable property data)  
Other Sources: DSL\*\*, EINECS\*\*, TSCA\*\*  
(\*\*Enter CHEMLIST File for up-to-date regulatory information)

Au

\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\*

138772 REFERENCES IN FILE CA (1907 TO DATE)  
4092 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA  
138913 REFERENCES IN FILE CAPLUS (1907 TO DATE)  
1 REFERENCES IN FILE CAOLD (PRIOR TO 1967)

L4 ANSWER 6 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN

RN 3458-28-4 REGISTRY

CN D-Mannose (9CI) (CA INDEX NAME)

OTHER CA INDEX NAMES:

CN Mannose, D- (8CI)

OTHER NAMES:

CN (+)-Mannose

CN Carubinose

CN D(+)-Mannose

CN Mannose

CN NSC 26247

CN Seminose

AR 530-26-7

FS STEREOSEARCH

DR 147-74-0

MF C6 H12 O6

CI COM

LC STN Files: ADISNEWS, AGRICOLA, BEILSTEIN\*, BIOBUSINESS, BIOSIS,  
BIOTECHNO, CA, CABA, CAOLD, CAPLUS, CASREACT, CEN, CHEMCATS,  
CHEMINFORMRX, CHEMLIST, CIN, CSCHM, DETHERM\*, EMBASE, GMELIN\*, HODOC\*,  
IFICDB, IFIPAT, IFIUDB, IPA, MEDLINE, MRCK\*, MSDS-OHS, NAPRALERT,  
NIOSHTIC, PIRA, PROMT, SPECINFO, SYNTHLINE, TOXCENTER, TULSA, USPAT2,  
USPATFULL


(\*File contains numerically searchable property data)

Other Sources: DSL\*\*, EINECS\*\*, TSCA\*\*

(\*\*Enter CHEMLIST File for up-to-date regulatory information)

Absolute stereochemistry. Rotation (+).





\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\*

13669 REFERENCES IN FILE CA (1907 TO DATE)  
 623 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA  
 13685 REFERENCES IN FILE CAPLUS (1907 TO DATE)  
 7 REFERENCES IN FILE CAOLD (PRIOR TO 1967)

L4 ANSWER 7 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN

RN 2438-80-4 REGISTRY

CN L-Galactose, 6-deoxy- (9CI) (CA INDEX NAME)

OTHER CA INDEX NAMES:

CN Fucose, L- (8CI)

OTHER NAMES:

CN (-)-Fucose

CN (-)-L-Fucose

CN 6-Deoxy-L-galactose

CN 6-Desoxygalactose

CN Fucose

CN L-(-)-Fucose

CN L-Fucose

CN L-Galactomethylose

AR 87-96-7, 3713-31-3

FS STEREOSEARCH

MF C6 H12 O5


CI COM

LC STN Files: AGRICOLA, ANABSTR, BEILSTEIN\*, BIOBUSINESS, BIOSIS,  
 BIOTECHNO, CA, CABA, CAPLUS, CASREACT, CEN, CHEMCATS, CHEMLIST, CIN,  
 CSCHEM, DDFU, DRUGU, EMBASE, MEDLINE, MRCK\*, MSDS-OHS, NAPRALERT,  
 NIOSHTIC, PIRA, PROMT, TOXCENTER, TULSA, USPAT2, USPATFULL  
 (\*File contains numerically searchable property data)

Other Sources: EINECS\*\*, NDSL\*\*, TSCA\*\*

(\*\*Enter CHEMLIST File for up-to-date regulatory information)

Absolute stereochemistry.



\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\*

5035 REFERENCES IN FILE CA (1907 TO DATE)  
240 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA  
5039 REFERENCES IN FILE CAPLUS (1907 TO DATE)

L4 ANSWER 8 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN  
RN 151-51-9 REGISTRY  
CN Methanediimine (9CI) (CA INDEX NAME)  
OTHER CA INDEX NAMES:  
CN Carbodiimide (6CI, 7CI, 8CI)  
OTHER NAMES:  
CN Stabilisator 9000  
FS 3D CONCORD  
MF C H2 N2  
CI COM  
LC STN Files: AGRICOLA, BEILSTEIN\*, BIOBUSINESS, BIOSIS, BIOTECHNO, CA,  
CAOLD, CAPLUS, CASREACT, CEN, CHEMCATS, CHEMLIST, CIN, CSNB, EMBASE,  
GMELIN\*, IFICDB, IFIPAT, IFIUDB, PIRA, PROMT, TOXCENTER, USPAT2,  
USPATFULL  
(\*File contains numerically searchable property data)

HN=C=NH


\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\*

733 REFERENCES IN FILE CA (1907 TO DATE)  
199 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA  
734 REFERENCES IN FILE CAPLUS (1907 TO DATE)  
4 REFERENCES IN FILE CAOLD (PRIOR TO 1967)

L4 ANSWER 9 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN  
RN 131-48-6 REGISTRY  
CN Neuraminic acid, N-acetyl- (9CI) (CA INDEX NAME)  
OTHER CA INDEX NAMES:  
CN D-glycero-D-galacto-2-Nonulosonic acid, 5-(acetylamino)-3,5-dideoxy-  
CN D-glycero-D-galacto-Nonulosonic acid, 5-acetamido-3,5-dideoxy- (8CI)  
CN Lactaminic acid (7CI)  
OTHER NAMES:  
CN 5-N-Acetyl-D-neuraminic acid  
CN 5-N-Acetylneuraminic acid  
CN Aceneuramic acid  
CN Acetylneuraminic acid  
CN N-Acetyl-D-neuraminic acid  
CN N-Acetylneuramic acid  
CN N-Acetylneuraminic acid  
CN N-Acetylsialic acid  
CN NANA  
FS STEREOSEARCH  
DR 6918-20-3, 11032-36-3, 14752-56-8, 5977-25-3, 6225-16-7  
MF C11 H19 N O9  
CI COM  
LC STN Files: ADISINSIGHT, AGRICOLA, ANABSTR, BEILSTEIN\*, BIOBUSINESS,  
BIOSIS, BIOTECHNO, CA, CANCERLIT, CAOLD, CAPLUS, CASREACT, CEN,  
CHEMCATS, CHEMLIST, CIN, CSCHM, DDFU, DRUGU, EMBASE, IFICDB, IFIPAT,  
IFIUDB, IPA, MEDLINE, MSDS-OHS, NIOSHTIC, PHAR, PROMT, PROUSDDR,

SPECINFO, SYNTHLINE, TOXCENTER, USAN, USPAT2, USPATFULL  
 (\*File contains numerically searchable property data)  
 Other Sources: EINECS\*\*, NDSL\*\*, TSCA\*\*  
 (\*\*Enter CHEMLIST File for up-to-date regulatory information)

Absolute stereochemistry.



\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\*

2516 REFERENCES IN FILE CA (1907 TO DATE)  
 161 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA  
 2518 REFERENCES IN FILE CAPLUS (1907 TO DATE)  
 1 REFERENCES IN FILE CAOLD (PRIOR TO 1967)

L4 ANSWER 10 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN

RN 59-23-4 REGISTRY

CN D-Galactose (9CI) (CA INDEX NAME)

OTHER CA INDEX NAMES:

CN Galactose, D- (8CI)

OTHER NAMES:

CN (+)-Galactose

CN D-(+)-Galactose

CN Galactose

FS STEREOSEARCH

DR 147-76-2, 3812-56-4, 400876-94-0

MF C6 H12 O6

CI COM


LC STN Files: ADISNEWS, AGRICOLA, ANABSTR, BEILSTEIN\*, BIOBUSINESS, BIOSIS, BIOTECHNO, CA, CABA, CAOLD, CAPLUS, CASREACT, CBNB, CEN, CHEMCATS, CHEMINFORMRX, CHEMLIST, CIN, CSCHEM, DDFU, DETHERM\*, DIOGENES, DRUGU, EMBASE, GMELIN\*, HODOC\*, IFICDB, IFIPAT, IFIUDB, IPA, MEDLINE, MRCK\*, MSDS-OHS, NAPRALERT, NIOSHTIC, PIRA, PROMT, RTECS\*, SPECINFO, SYNTHLINE, TOXCENTER, TULSA, USPAT2, USPATFULL, VETU

(\*File contains numerically searchable property data)

Other Sources: DSL\*\*, EINECS\*\*, TSCA\*\*

(\*\*Enter CHEMLIST File for up-to-date regulatory information)

Absolute stereochemistry. Rotation (+).



**\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\***

21536 REFERENCES IN FILE CA (1907 TO DATE)  
761 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA  
21559 REFERENCES IN FILE CAPLUS (1907 TO DATE)  
2 REFERENCES IN FILE CAOLD (PRIOR TO 1967)

L4 ANSWER 11 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN

RN 58-86-6 REGISTRY

CN D-Xylose (9CI) (CA INDEX NAME)

OTHER CA INDEX NAMES:

CN Xylose, D- (8CI)

OTHER NAMES:

CN (+)-Xylose

CN D-(+)-Xylose

CN Wood sugar

CN Xylose

FS STEREOSEARCH

DR 133-56-2, 141492-19-5

MF C5 H10 O5

CI COM


LC STN Files: ADISNEWS, AGRICOLA, ANABSTR, BEILSTEIN\*, BIOBUSINESS, BIOSIS, BIOTECHNO, CA, CABA, CAOLD, CAPLUS, CASREACT, CBNB, CHEMCATS, CHEMINFORMRX, CHEMLIST, CIN, CSCHM, DDFU, DETHERM\*, DIOGENES, DRUGU, EMBASE, GMELIN\*, HODOC\*, HSDB\*, IFICDB, IFIPAT, IFIUDB, IPA, MEDLINE, MRCK\*, MSDS-OHS, NAPRALERT, NIOSHTIC, PIRA, PROMT, RTECS\*, SPECINFO, SYNTHLINE, TOXCENTER, TULSA, USAN, USPAT2, USPATFULL, VETU, VTB

(\*File contains numerically searchable property data)

Other Sources: DSL\*\*, EINECS\*\*, TSCA\*\*

(\*\*Enter CHEMLIST File for up-to-date regulatory information)

Absolute stereochemistry.

**\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\***

13871 REFERENCES IN FILE CA (1907 TO DATE)  
317 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA  
13882 REFERENCES IN FILE CAPLUS (1907 TO DATE)  
5 REFERENCES IN FILE CAOLD (PRIOR TO 1967)

L4 ANSWER 12 OF 12 REGISTRY COPYRIGHT 2004 ACS on STN

RN 50-99-7 REGISTRY

CN D-Glucose (8CI, 9CI) (CA INDEX NAME)


OTHER NAMES:

CN (+)-Glucose

CN Anhydrous dextrose

CN Cartose  
 CN Cerelose  
 CN Cerelose 2001  
 CN Clearsweet 95  
 CN Clintose L  
 CN Corn sugar  
 CN CPC hydrate  
 CN D(+)-Glucose  
 CN D-glucose  
 CN Dextropur  
 CN Dextrose  
 CN Dextrosol  
 CN Glucodin  
 CN Glucolin  
 CN Glucose  
 CN Glucosteril  
 CN Goldsugar  
 CN Grape sugar  
 CN Maxim Energy Gel  
 CN Roferose ST  
 CN Staleydex 111  
 CN Staleydex 130  
 CN Staleydex 333  
 CN Sugar, grape  
 CN Tabfine 097(HS)  
 CN Vadex  
 FS STEREOSEARCH  
 DR 8012-24-6, 8030-23-7, 162222-91-5, 165659-51-8, 50933-92-1, 80206-31-1  
 MF C6 H12 O6  
 CI COM  
 LC STN Files: ADISNEWS, AGRICOLA, ANABSTR, AQUIRE, BEILSTEIN\*, BIOBUSINESS,  
 BIOSIS, BIOTECHNO, CA, CABA, CANCERLIT, CAOLD, CAPLUS, CASREACT, CBNB,  
 CEN, CHEMCATS, CHEMINFORMRX, CHEMLIST, CHEMSAFE, CIN, CSCHEM, CSNB,  
 DDFU, DETHERM\*, DIOGENES, DIPPR\*, DRUGU, EMBASE, GMELIN\*, HSDB\*, IFICDB,  
 IFIPAT, IFIUDB, IMSCOSEARCH, IPA, MEDLINE, MRCK\*, MSDS-OHS, NAPRALERT,  
 NIOSHTIC, PDLCOM\*, PIRA, PROMT, PS, RTECS\*, SPECINFO, TOXCENTER, TULSA,  
 ULIDAT, USAN, USPAT2, USPATFULL, VETU, VTB  
 (\*File contains numerically searchable property data)  
 Other Sources: DSL\*\*, EINECS\*\*, TSCA\*\*  
 (\*\*Enter CHEMLIST File for up-to-date regulatory information)

Absolute stereochemistry.




\*\*PROPERTY DATA AVAILABLE IN THE 'PROP' FORMAT\*\*

167083 REFERENCES IN FILE CA (1907 TO DATE)  
 2199 REFERENCES TO NON-SPECIFIC DERIVATIVES IN FILE CA  
 167254 REFERENCES IN FILE CAPLUS (1907 TO DATE)

14 REFERENCES IN FILE CAOLD (PRIOR TO 1967)

L5 STR



## NODE ATTRIBUTES:

DEFAULT MLEVEL IS ATOM

DEFAULT ECLEVEL IS LIMITED

## GRAPH ATTRIBUTES:


RING(S) ARE ISOLATED OR EMBEDDED

NUMBER OF NODES IS 6

## STEREO ATTRIBUTES: NONE

L8 9777 SEA FILE=REGISTRY SSS FUL L5

L9 STR



-OEESEECONHNH-

## NODE ATTRIBUTES:

DEFAULT MLEVEL IS ATOM

DEFAULT ECLEVEL IS LIMITED

## GRAPH ATTRIBUTES:

RING(S) ARE ISOLATED OR EMBEDDED

NUMBER OF NODES IS 10

## STEREO ATTRIBUTES: NONE

L10 17 SEA FILE=REGISTRY SUB=L8 SSS FUL L9

L18 10 SEA FILE=HCAPLUS ABB=ON PLU=ON L10

L19 2 SEA FILE=HCAPLUS ABB=ON PLU=ON L18 AND (IMMOBIL? OR COAT? OR ATTACH? OR BIOSENS? OR BIOCHIP? OR BIO?(2A)?SENS?)

L20 10 SEA FILE=HCAPLUS ABB=ON PLU=ON L18 OR L19

=&gt; d l20 ibib abs hitind hitstr 1-10

L20 ANSWER 1 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1996:711448 HCAPLUS

DOCUMENT NUMBER: 126:89666

TITLE: Synthesis of ligands related to the Vibrio cholerae O-specific antigen. Part 12. Synthesis of eight glycosides of hexasaccharide fragments representing the terminus of the O-polysaccharide of Vibrio cholerae O:1, serotype Inaba and Ogawa, bearing aglycons suitable for linking to proteins

AUTHOR(S): Ogawa, Yuji; Lei, Ping-sheng; Kovac, Pavol

CORPORATE SOURCE: Natl. Inst. Health, NIDDK, Bethesda, MD, 20892-0815, USA

SOURCE: Carbohydrate Research (1996) 293(2), 173-194

CODEN: CRBRAT; ISSN: 0008-6215

PUBLISHER: Elsevier  
DOCUMENT TYPE: Journal  
LANGUAGE: English

AB The title substances were prepared from intermediate, fully acetylated  $\alpha$ -trimethylsilylethyl (SE) glycosides. The latter were assembled in a blockwise manner, using as the glycosyl donor the  $\alpha$ -glycosyl chloride of a disaccharide bearing two 4-azido-4-deoxy functions. Next, the azido groups in the assembled hexasaccharides were converted to the corresponding amines, and these were acylated with 4-O-benzyl-3-deoxy-L-glycero-tetronic acid in the presence of a water-soluble carbodiimide. The Se glycoside were then transformed to glycosyl imidates, and these were coupled with Me 6-hydroxyhexanoate or Me 2-(2-hydroxyethylthio)propionate. The aglycons in the glycosides thus obtained were then converted to the corresponding carboxylic acids or acyl hydrazides. Such compds. are suitable for linking to proteins to obtain neoglycoproteins.

CC 33-7 (Carbohydrates)

IT 4547-43-7P, Methyl 6-hydroxyhexanoate 185248-27-5P 185248-28-6P  
185248-29-7P 185248-30-0P 185248-54-8P 185248-55-9P  
185248-56-0P 185248-57-1P 185248-58-2P

RL: SPN (Synthetic preparation); PREP (Preparation)  
(synthesis of eight glycosides of hexasaccharide fragments representing the terminus of the O-polysaccharide of *Vibrio cholerae*)

IT 185248-54-8P 185248-58-2P

RL: SPN (Synthetic preparation); PREP (Preparation)  
(synthesis of eight glycosides of hexasaccharide fragments representing the terminus of the O-polysaccharide of *Vibrio cholerae*)


RN 185248-54-8 HCAPLUS

CN Propanoic acid, 3-[[2-[[O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl-(1 $\rightarrow$ 2)-O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl-(1 $\rightarrow$ 2)-O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl-(1 $\rightarrow$ 2)-O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl-(1 $\rightarrow$ 2)-O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl]oxy]ethylthio]-, hydrazide (9CI) (CA INDEX NAME)


Absolute stereochemistry. Rotation (-).




PAGE 1-A



PAGE 1-B



PAGE 2-A




PAGE 2-B




RN 185248-58-2 HCAPLUS  
 CN Propanoic acid, 3-[[2-[[[O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]-2-O-methyl- $\alpha$ -D-mannopyranosyl-(1 $\rightarrow$ 2)-O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl-(1 $\rightarrow$ 2)-O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl-(1 $\rightarrow$ 2)-O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl-(1 $\rightarrow$ 2)-O-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl-(1 $\rightarrow$ 2)-4,6-dideoxy-4-[[[(2S)-2,4-dihydroxy-1-oxobutyl]amino]- $\alpha$ -D-mannopyranosyl]oxy]ethyl]thio]-, hydrazide (9CI) (CA INDEX NAME)

Absolute stereochemistry. Rotation (-).


PAGE 1-A




PAGE 1-B



PAGE 2-A



PAGE 2-B



REFERENCE COUNT: 24 THERE ARE 24 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L20 ANSWER 2 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1996:504622 HCAPLUS

DOCUMENT NUMBER: 125:276374

TITLE: Synthesis of ligands related to *Vibrio cholerae* O-specific antigen. II. Synthesis of four glycosides of a disaccharide fragment representing the terminus of the O-polysaccharide of *Vibrio cholerae* O:1, serotype Inaba, bearing aglycons suitable for linking to proteins

AUTHOR(S): Ogawa, Yuji; Lei, Ping-sheng; Kovac, Pavol

CORPORATE SOURCE: NIDDK, Nat. Inst. Health, Bethesda, MD, 20892-0815, USA

SOURCE: Carbohydrate Research (1996), 288, 85-98  
CODEN: CRBRAT; ISSN: 0008-6215

PUBLISHER: Elsevier

DOCUMENT TYPE: Journal

LANGUAGE: English

AB Me 4-azido-3-O-benzyl-4,6-dideoxy- $\alpha$ -D-mannopyranoside was converted into disaccharide glycosides fragment of the terminus of the O-polysaccharide of *Vibrio cholerae* O:1.

CC 33-7 (Carbohydrates)

IT 182273-70-7P 182273-73-0P

RL: SPN (Synthetic preparation); PREP (Preparation)  
(preparation of glycosides of a disaccharide fragment representing the terminus of O-polysaccharide of *Vibrio cholerae* O:1)


IT 182273-73-0P

RL: SPN (Synthetic preparation); PREP (Preparation)  
(preparation of glycosides of a disaccharide fragment representing the

terminus of O-polysaccharide of *Vibrio cholerae* O:1)

RN 182273-73-0 HCAPLUS  
 CN Propanoic acid, 3-[[[2-[[[4,6-dideoxy-2-O-[4,6-dideoxy-4-[(2,4-dihydroxy-1-oxobutyl)amino]- $\alpha$ -D-mannopyranosyl]-4-[(2,4-dihydroxy-1-oxobutyl)amino]- $\alpha$ -D-mannopyranosyl]oxy]ethyl]thio]-, hydrazide, [2(S),4(S)]- (9CI) (CA INDEX NAME)

Absolute stereochemistry. Rotation (+).



L20 ANSWER 3 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN  
 ACCESSION NUMBER: 1995:241713 HCAPLUS  
 DOCUMENT NUMBER: 122:75648  
 TITLE: A carbohydrate **biosensor** surface for the detection of uropathogenic bacteria  
 AUTHOR(S): Nilsson, Kurt G. I.; Mandenius, Carl-Fredrik  
 CORPORATE SOURCE: Lund, S-223 70, Swed.  
 SOURCE: Bio/Technology (1994), 12(13), 1376-8  
 CODEN: BTCHDA; ISSN: 0733-222X  
 PUBLISHER: Nature Publishing Co.  
 DOCUMENT TYPE: Journal  
 LANGUAGE: English

AB We have developed a new surface for use in **biosensors** that is based on a gold plate covered with a specific carbohydrate receptor structure. The carbohydrate, Gal $\alpha$ 1-4Gal, was bound covalently via a thioalkylcarboxy-spacer, or adsorbed as a neoglycoprotein, to a two-dimensional gold surface. Both types of surfaces showed high specificity in the binding of the uropathogenic bacteria P-fimbriated *Escherichia coli* compared to the binding of non-infectious bacteria. The signal to noise ratio is sufficiently high to allow specific detection of the bacteria in **biosensor** applications.

CC 9-1 (Biochemical Methods)

ST carbohydrate **biosensor** surface uropathogenic bacteria detection

IT *Escherichia coli*

(P-fimbriated; carbohydrate **biosensor** surface for the detection of uropathogenic bacteria)

IT Albumins, uses


RL: DEV (Device component use); PEP (Physical, engineering or chemical process); PROC (Process); USES (Uses)

(bovine, carbohydrate **biosensor** surface for the detection of uropathogenic bacteria)

IT **Biosensors**


- (carbohydrate **biosensor** surface for the detection of uropathogenic bacteria)
- IT Carbohydrates and Sugars, uses  
 RL: DEV (Device component use); USES (Uses)  
 (carbohydrate **biosensor** surface for the detection of uropathogenic bacteria)
- IT Bacteria  
 (uropathogenic, carbohydrate **biosensor** surface for the detection of uropathogenic bacteria)
- IT 7440-57-5, Gold, uses  
 RL: DEV (Device component use); USES (Uses)  
 (carbohydrate **biosensor** surface for the detection of uropathogenic bacteria)
- IT 160294-57-5 160294-57-5D, conjugate with bovine serum albumin  
 RL: DEV (Device component use); PEP (Physical, engineering or chemical process); PROC (Process); USES (Uses)  
 (carbohydrate **biosensor** surface for the detection of uropathogenic bacteria)
- IT 160294-57-5 160294-57-5D, conjugate with bovine serum albumin  
 RL: DEV (Device component use); PEP (Physical, engineering or chemical process); PROC (Process); USES (Uses)  
 (carbohydrate **biosensor** surface for the detection of uropathogenic bacteria)
- RN 160294-57-5 HCAPLUS  
 CN Propanoic acid, 3-[[2-[(4-O- $\alpha$ -D-galactopyranosyl- $\beta$ -D-galactopyranosyl)oxy]ethyl]thio]-, hydrazide (9CI) (CA INDEX NAME)

Absolute stereochemistry.



- RN 160294-57-5 HCAPLUS  
 CN Propanoic acid, 3-[[2-[(4-O- $\alpha$ -D-galactopyranosyl- $\beta$ -D-galactopyranosyl)oxy]ethyl]thio]-, hydrazide (9CI) (CA INDEX NAME)

Absolute stereochemistry.



L20 ANSWER 4 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1989:587593 HCAPLUS

DOCUMENT NUMBER: 111:187593

TITLE: Amine derivatives of anthracycline antibiotics and antibody conjugates thereof and their preparation and use in treatment of cellular disorders

INVENTOR(S): King, Dalton H.; Coughlin, Daniel J.; Rodwell, John Dennis; Lopes, Anthony Dwight; Radcliffe, Robert David

PATENT ASSIGNEE(S): Cytogen Corp., USA

SOURCE: Eur. Pat. Appl., 41 pp.

CODEN: EPXXDW

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 5

PATENT INFORMATION:


PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
EP 294294	A2	19881207	EP 1988-401353	19880603
EP 294294	A3	19900530		
EP 294294	B1	19950517		
R: AT, BE, CH, DE, ES, FR, GB, GR, IT, LI, LU, NL, SE				
US 4950738	A	19900821	US 1987-58440	19870605
US 5162512	A	19921110	US 1988-199549	19880527
WO 8809823	A1	19881215	WO 1988-US1909	19880603
W: AU, DK, FI, JP				
AU 8819549	A1	19890104	AU 1988-19549	19880603
ZA 8803956	A	19890222	ZA 1988-3956	19880603
JP 02500749	T2	19900315	JP 1988-505203	19880603
ES 2074055	T3	19950901	ES 1988-401353	19880603
DK 8900511	A	19890203	DK 1989-511	19890203
FI 8900559	A	19890206	FI 1989-559	19890206
PRIORITY APPLN. INFO.:				
			US 1987-58440	19870605
			US 1988-199549	19880527
			US 1982-356315	19820309
			US 1984-650375	19840913
			US 1984-650754	19840913
			WO 1988-US1909	19880603

OTHER SOURCE(S): MARPAT 111:187593

AB Amine derivs. of antineoplastic anthracycline antibiotics (e.g. hydrazine, hydrazide, phenylhydrazine, etc. derivs. of daunorubicin, doxorubicin,

carminocyan, etc.) are prepared and covalently **attached** to an antibody or antibody fragment for treatment of cellular disorders, especially neoplasms. Adriamycin-HCl was reacted with adipic dihydrazide and the product (ADR-ADH) was conjugated with a murine monoclonal antibody B72.3 specific for an antigen of human adenocarcinoma (the oligosaccharide moiety of the antibody had been oxidized with NaIO<sub>4</sub>). Tumor (human colon adenocarcinoma, BL/CX-3) growth in nude mice treated i.v. with 6 µg ADR-ADH-B72.3 conjugate was significantly inhibited compared to the untreated group. The tumor inhibitory effect was equivalent to that seen in animals receiving 200 µg ADR alone. The inhibition lasted beyond the end of the treatment.

IC ICM C07H015-252  
ICS A61K031-70; C07K015-00; A61K039-395; A61K047-00  
CC 1-6 (Pharmacology)  
Section cross-reference(s): 9, 15  
IT 23214-92-8DP, pentaglutamylhydrazide derivs., antibody conjugates  
123105-65-7DP, antibody conjugates 123105-66-8DP, antibody conjugates  
123105-67-9DP, antibody conjugates 123105-68-0DP, antibody conjugates  
123105-69-1DP, antibody conjugates 123105-70-4DP, antibody conjugates  
123105-71-5DP, antibody conjugates 123105-72-6DP, antibody conjugates  
123105-73-7DP, antibody conjugates 123105-74-8DP, antibody conjugates  
123105-75-9DP, antibody conjugates 123106-26-3DP, antibody conjugates  
**123129-58-8DP, antibody conjugates**  
RL: SPN (Synthetic preparation); PREP (Preparation)  
(preparation of, as neoplasm inhibitors)  
IT **123129-58-8DP, antibody conjugates**  
RL: SPN (Synthetic preparation); PREP (Preparation)  
(preparation of, as neoplasm inhibitors)  
RN 123129-58-8 HCAPLUS  
CN Propanoic acid, 3-[[2-[4-[(3-amino-2,3,6-trideoxy-α-L-lyxo-hexopyranosyl)oxy]-1,2,3,4,6,11-hexahydro-2,5,12-trihydroxy-7-methoxy-6,11-dioxo-2-naphthacenyl]-2-oxoethyl]thio]-, hydrazide, (2S-cis)- (9CI) (CA INDEX NAME)



L20 ANSWER 5 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1986:573419 HCAPLUS

DOCUMENT NUMBER: 105:173419

TITLE: Polyurethane ionic polymers with sulfide and semicarbazide groups in the macrochain



AUTHOR(S): Sukhorukova, S. A.; Navrotskaya, R. P.; Grekov, A. P.  
 CORPORATE SOURCE: Inst. Khim. Vysokomol. Soedin., Kiev, USSR  
 SOURCE: Ukrainskii Khimicheskii Zhurnal (Russian Edition)  
 (1986), 52(5), 540-3  
 CODEN: UKZHAU; ISSN: 0041-6045

DOCUMENT TYPE: Journal  
 LANGUAGE: Russian

AB S-containing polyurethane ionenes prepared from polytetramethylene glycol, TDI or MDI, (H<sub>2</sub>NNHCO)<sub>2</sub>Z [Z= CH<sub>2</sub>CH(SR), O(CH<sub>2</sub>CH<sub>2</sub>SCH<sub>2</sub>CH<sub>2</sub>)<sub>2</sub>; R = alkyl], and (HOCH<sub>2</sub>CH<sub>2</sub>)<sub>2</sub>NMe with subsequent quaternization with HCl and optional alkylation with EtBr or Me<sub>2</sub>SO<sub>4</sub> were more UV resistant than their nonsulfur analogs. No changes were observed in tensile strength, elongation, or viscosity after 100-h UV exposure. The photostabilizing influence of S was observed also for the nonquaternized polyurethanes. The S-containing fragments apparently generated weakly active radicals of the type RS·, which were capable of terminating the kinetic chain of degradation

CC 36-5 (Physical Properties of Synthetic High Polymers)

IT 72186-71-1D, quaternized 72186-72-2 72186-72-2D, quaternized  
 72186-74-4 72186-74-4D, quaternized 72196-93-1 72196-93-1D,  
 quaternized 104935-10-6D, quaternized 104935-11-7D,  
 quaternized 104985-03-7 104985-04-8 104985-05-9 104985-06-0  
 104985-08-2 104985-09-3 104985-10-6  
 RL: PRP (Properties)  
 (UV stability of, physicomech. properties in relation to)


IT 104935-11-7D, quaternized 104985-08-2  
 104985-09-3 104985-10-6  
 RL: PRP (Properties)  
 (UV stability of, physicomech. properties in relation to)

RN 104935-11-7 HCAPLUS

CN Propanoic acid, 3,3'-[oxybis(2,1-ethanediylthio)]bis-, dihydrazide,  
 polymer with α-hydro-ω-hydroxypoly(oxy-1,4-butanediyl) and  
 2,2'-(methyylimino)bis[ethanol] (9CI) (CA INDEX NAME)


CM 1

CRN 92268-36-5  
 CMF C10 H22 N4 O3 S2




CM 2

CRN 25190-06-1  
 CMF (C<sub>4</sub> H<sub>8</sub> O)<sub>n</sub> H<sub>2</sub> O  
 CCI PMS



CM 3

CRN 105-59-9  
CMF C5 H13 N O2




RN 104985-08-2 HCAPLUS  
CN Propanoic acid, 3,3'-[oxybis(2,1-ethanediylthio)]bis-, dihydrazide, polymer with 1,3-diisocyanatomethylbenzene,  $\alpha$ -hydro- $\omega$ -hydroxypoly(oxy-1,4-butanediyl) and 2,2'-(methyylimino)bis[ethanol], hydrochloride (9CI) (CA INDEX NAME)

CM 1

CRN 104985-07-1  
CMF (C10 H22 N4 O3 S2 . C9 H6 N2 O2 . C5 H13 N O2 . (C4 H8 O)n H2 O)x  
CCI PMS


CM 2

CRN 92268-36-5  
CMF C10 H22 N4 O3 S2



CM 3


CRN 26471-62-5  
CMF C9 H6 N2 O2  
CCI IDS



D1-Me

CM 4


CRN 25190-06-1  
CMF (C4 H8 O)n H2 O  
CCI PMS



CM 5

CRN 105-59-9

CMF C5 H13 N O2




RN 104985-09-3 HCAPLUS

CN Propanoic acid, 3,3'-[oxybis(2,1-ethanedithio)]bis-, dihydrazide, polymer with 1,3-diisocyanatomethylbenzene,  $\alpha$ -hydro- $\omega$ -hydroxypoly(oxy-1,4-butanediyl) and 2,2'-(methylimino)bis[ethanol], compd. with bromoethane, hydrochloride (9CI) (CA INDEX NAME)

CM 1

CRN 74-96-4

CMF C2 H5 Br



CM 2

CRN 104985-07-1


CMF (C10 H22 N4 O3 S2 . C9 H6 N2 O2 . C5 H13 N O2 . (C4 H8 O)n H2 O)x

CCI PMS

CM 3

CRN 92268-36-5

CMF C10 H22 N4 O3 S2




CM 4

CRN 26471-62-5

CMF C9 H6 N2 O2

CCI IDS




D1-- Me

CM 5

CRN 25190-06-1

CMF (C4 H8 O)n H2 O


CCI PMS



CM 6

CRN 105-59-9

CMF C5 H13 N O2




RN 104985-10-6 HCAPLUS

CN Propanoic acid, 3,3'-[oxybis(2,1-ethanediylthio)]bis-, dihydrazide, polymer with 1,3-diisocyanatomethylbenzene,  $\alpha$ -hydro- $\omega$ -hydroxypoly(oxy-1,4-butanediyl) and 2,2'-(methyylimino)bis[ethanol], compd. with dimethyl sulfate, hydrochloride (9CI) (CA INDEX NAME)

CM 1

CRN 77-78-1

CMF C2 H6 O4 S



CM 2

CRN 104985-07-1


CMF (C10 H22 N4 O3 S2 . C9 H6 N2 O2 . C5 H13 N O2 . (C4 H8 O)n H2 O)x

CCI PMS

CM 3

CRN 92268-36-5

CMF C10 H22 N4 O3 S2




CM 4

CRN 26471-62-5

CMF C9 H6 N2 O2

CCI IDS




D1-Me

CM 5

CRN 25190-06-1

CMF (C4 H8 O)<sub>n</sub> H2 O


CCI PMS



CM 6


CRN 105-59-9

CMF C5 H13 N O2




L20 ANSWER 6 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN  
 ACCESSION NUMBER: 1986:169828 HCAPLUS

DOCUMENT NUMBER: 104:169828  
 TITLE: Effect of UV irradiation on sulfur-containing poly(urethane semicarbazide)s  
 AUTHOR(S): Sukhorukova, S. A.; Navrotskaya, R. P.; Grekov, A. P.; Fedorenko, O. M.  
 CORPORATE SOURCE: Inst. Khim. Vysokomol. Soedin., Moscow, USSR  
 SOURCE: Vysokomolekulyarnye Soedineniya, Seriya A (1986), 28(1), 111-16  
 CODEN: VYSAAF; ISSN: 0507-5475  
 DOCUMENT TYPE: Journal  
 LANGUAGE: Russian  
 AB S-containing polyurethane-polysemicarbazide (PUS) rubbers prepared from polytetramethylene glycol, MDI, and thiodicarboxylic acid dihydrazides, containing a sulfide group in the main chain or in a side group were more stable to UV degradation than similar polyurethanes not containing semicarbazide or S groups. The stabilizing influence of the semicarbazide groups was attributed to its antioxidant activity resulting from the presence of replaceable H atoms. The S-containing fragments stabilized through formation of RS radicals capable of terminating the kinetic chain of degradation. Photodegrdn. of the S-containing PUS occurred only in the first 5-10 h of irradiation, after which the mech. properties increased in value then became constant with further irradiation  
 CC 39-7 (Synthetic Elastomers and Natural Rubber)  
 IT 9018-04-6 52484-70-5 95410-92-7 95410-93-8 101909-11-9  
 101909-12-0 101909-13-1 **101909-14-2**  
 RL: PRP (Properties)  
 (stability of, to photodegrdn. by UV irradiation, structure in relation to)  
 IT **101909-14-2**  
 RL: PRP (Properties)  
 (stability of, to photodegrdn. by UV irradiation, structure in relation to)  
 RN 101909-14-2 HCAPLUS  
 CN Propanoic acid, 3,3'-[oxybis(2,1-ethanediylthio)]bis-, dihydrazide, polymer with  $\alpha$ -hydro- $\omega$ -hydroxypoly(oxy-1,4-butanediyl) and 1,1'-methylenebis[4-isocyanatobenzene] (9CI) (CA INDEX NAME)  
 CM 1  
 CRN 92268-36-5  
 CMF C10 H22 N4 O3 S2




CM 2

CRN 25190-06-1  
 CMF (C4 H8 O)<sub>n</sub> H2 O  
 CCI PMS



CM 3

CRN 101-68-8  
CMF C15 H10 N2 O2



L20 ANSWER 7 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1985:168202 HCAPLUS

DOCUMENT NUMBER: 102:168202

TITLE: Flexible polyurethanes

INVENTOR(S): Sukhorukova, S. A.; Navrotskaya, R. P.; Grekov, A. P.;  
Tanchuk, Yu. V.; Kornienko, A. A.

PATENT ASSIGNEE(S): Institute of the Chemistry of High-Molecular-Weight  
Compounds, Academy of Sciences, Ukrainian S.S.R., USSR

SOURCE: U.S.S.R. From: Otkrytiya, Izobret. 1984, (48), 88-9.

CODEN: URXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Russian

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
SU 1131886	A1	19841230	SU 1983-3587911	19830309

PRIORITY APPLN. INFO.: SU 1983-3587911 19830309

AB Flexible polyurethanes are prepared by reacting oligoesters and diisocyanates and then adding chain extenders at a mol. ratio of oligomers/diisocyanates/chain extenders of 1:2-2.2:0.9-1.2. The UV light and low-temperature resistance of the polyurethanes is increased by using as a chain extender 1,5-diethylene oxide-S,S'-bis(mercaptopropionic acid) dihydrazide have the formula  $H_2NHNCOCH_2CH_2-S-CH_2CH_2-O-CH_2CH_2-S-CH_2CH_2CONHNH_2$ .

IC ICM C08G018-38


CC 39-12 (Synthetic Elastomers and Natural Rubber)

IT **92268-36-5DP**, polymers with polyesters and diisocyanates  
RL: PREP (Preparation)  
(oligomeric, rubber, manufacture of, light-and low-temperature-resistant)

IT **92268-36-5DP**, polymers with polyesters and diisocyanates  
RL: PREP (Preparation)  
(oligomeric, rubber, manufacture of, light-and low-temperature-resistant)

RN 92268-36-5 HCAPLUS

CN Propanoic acid, 3,3'-[oxybis(2,1-ethanediylthio)]bis-, dihydrazide (9CI)  
(CA INDEX NAME)



L20 ANSWER 8 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1984:552748 HCAPLUS

DOCUMENT NUMBER: 101:152748

TITLE: Thermal stability of sulfur-containing polyurethane semicarbazides

AUTHOR(S): Grekov, A. P.; Navrotskaya, R. P.; Zapunnaya, K. V.; Sukhorukova, S. A.; Fedorenko, O. M.

CORPORATE SOURCE: Inst. Khim. Vysokomol. Soedin., Kiev, USSR

SOURCE: Ukrainskii Khimicheskii Zhurnal (Russian Edition) (1984), 50(6), 659-63

CODEN: UKZHAU; ISSN: 0041-6045

DOCUMENT TYPE: Journal

LANGUAGE: Russian

AB Thermogravimetric anal., DTA, and IR spectroscopy was used to study the thermal properties of polyurethane semicarbazides (PUS) containing sulfide groups in the main or side chains, and the effect of sulfide fragments of dicarboxylic dihydrazides on the thermal stability (TS) of PUS. The PUS were prepared from oligomeric polytetramethylene glycol, TDI, and different dihydrazides in DMF solns. The oxidative TS of PUS was significantly higher than that of polyurethanes. The degradation activation energy (Ea) PUS was 106-325 kJ/mol of 240-510°. The highest Ea (325.1 kJ/mol) was observed for PUS with SC14H29 groups of the succinic dihydrazide. The TS of PUS depends on the chemical structure of the hydrazide, and especially on the structure of S-containing fragments.

CC 37-5 (Plastics Manufacture and Processing)

IT 90967-18-3 92268-33-2 92268-34-3 92268-35-4 92268-37-6

RL: USES (Uses)

(thermal and oxidative thermal stability of)

IT 92268-37-6

RL: USES (Uses)

(thermal and oxidative thermal stability of)


RN 92268-37-6 HCAPLUS

CN Propanoic acid, 3,3'-[oxybis(2,1-ethanedithio)]bis-, dihydrazide, polymer with 1,3-diisocyanatomethylbenzene and  $\alpha$ -hydro- $\omega$ -hydroxypoly(oxy-1,4-butanediyl) (9CI) (CA INDEX NAME)

CM 1

CRN 92268-36-5

CMF C10 H22 N4 O3 S2




CM 2

CRN 26471-62-5

CMF C9 H6 N2 O2

CCI IDS






D1-Me

CM 3

CRN 25190-06-1

CMF (C4 H8 O)<sub>n</sub> H2 O

CCI PMS



L20 ANSWER 9 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1984:423882 HCAPLUS

DOCUMENT NUMBER: 101:23882

TITLE: Glycosides and glycoconjugates

INVENTOR(S): Dahmen, Jan; Frejd, Torbjoern; Magnusson, Goeran; Noori, Ghazi

PATENT ASSIGNEE(S): Svenska Sockerfabriks AB, Swed.

SOURCE: Eur. Pat. Appl., 112 pp.

CODEN: EPXXDW

DOCUMENT TYPE: Patent


LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
EP 98252	A2	19840111	EP 1983-850176	19830621
EP 98252	A3	19840404		
EP 98252	B1	19890607		
R: AT, BE, CH, DE, FR, GB, IT, LI, NL, SE				
FI 8302254	A	19831224	FI 1983-2254	19830620
FI 78707	B	19890531		
FI 78707	C	19890911		
DK 8302842	A	19831224	DK 1983-2842	19830620
AT 43847	E	19890615	AT 1983-850176	19830621
JP 59025399	A2	19840209	JP 1983-111989	19830623
US 4675392	A	19870623	US 1984-673796	19841121
PRIORITY APPLN. INFO.:			SE 1982-3925	19820623
			US 1983-504154	19830614
			EP 1983-850176	19830621

GI



AB Glycosides of the general formula (sugar) $n$ O(CH $_2$ ) $m$ SRR1 [ $n$  = 1-10;  $m$  = 2-20; R = alkylene of  $\leq 25$  C atoms, arylene; R1 = H, CHO, NO $_2$ , NH $_2$ , OH, SH, CO $_2$ H, CO $_2$ Me, CO $_2$ Et, CONHNH $_2$ , CON $_3$ , CH(OR $_2$ ) $_2$ , (R $_2$  = C1-4 alkyl)] were prepared and some of them were converted into glycoconjugates. Thus, glycoside I (R $_3$  = Br), prepared from HO(CH $_2$ ) $_2$ Br and the corresponding sugar peracetate, was treated with HS(CH $_2$ ) $_2$ CO $_2$ Me to give I [R $_3$  = S(CH $_2$ ) $_2$ CO $_2$ Me]. The latter was deacetylated to give glycoside II (R $_4$  = OMe), which was treated with H $_2$ NNH $_2$  to give II (R $_4$  = NHNH $_2$ ), which was coupled to bovine serum albumin by the acyl azide method to give the corresponding conjugate II (R $_4$  = bovine serum albumin).

IC C07H015-04; C07H003-06; C07H003-04; C07G007-00; C07G017-00  
CC 33-3 (Carbohydrates)

Section cross-reference(s): 6

IT 90214-63-4P 90214-66-7P 90214-86-1P

RL: RCT (Reactant); SPN (Synthetic preparation); PREP (Preparation); RACT (Reactant or reagent)

(preparation and coupling of, with albumin)

IT 90214-63-4P 90214-66-7P 90214-86-1P


RL: RCT (Reactant); SPN (Synthetic preparation); PREP (Preparation); RACT (Reactant or reagent)

(preparation and coupling of, with albumin)

RN 90214-63-4 HCAPLUS

CN Propanoic acid, 3-[[2-( $\beta$ -D-galactopyranosyloxy)ethyl]thio]-, hydrazide (9CI) (CA INDEX NAME)


Absolute stereochemistry.



RN 90214-66-7 HCAPLUS

CN Propanoic acid, 3-[[2-[[2-O-(6-deoxy- $\alpha$ -L-galactopyranosyl)- $\beta$ -D-galactopyranosyl]oxy]ethyl]thio]-, hydrazide (9CI) (CA INDEX NAME)

Absolute stereochemistry.




RN 90214-86-1 HCAPLUS

CN Propanoic acid, 3-[[2-[(O-α-D-galactopyranosyl-(1→4)-O-β-D-galactopyranosyl)-(1→4)-β-D-glucopyranosyl]oxy]ethyl]thio]-, hydrazide (9CI) (CA INDEX NAME)

Absolute stereochemistry.

PAGE 1-A



PAGE 1-B

—NH<sub>2</sub>

L20 ANSWER 10 OF 10 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1981:41138 HCAPLUS

DOCUMENT NUMBER: 94:41138

TITLE: Studies on structure-activity relation of TAPHA-type compounds as monoamine oxidase inhibitors

AUTHOR(S): Wang, Yu-Ee; Xu, Fu-Ben; Chen, Chi-Hao; Jin, Guo-Zhang

CORPORATE SOURCE: Shanghai Inst. Mater. Med., Acad. Sin., Shanghai, Peop. Rep. China

SOURCE: Yaoxue Xuebao (1980), 15(3), 147-52

CODEN: YHHPAL; ISSN: 0513-4870

DOCUMENT TYPE: Journal

LANGUAGE: Chinese

AB The title compds. X[(CH<sub>2</sub>)<sub>n</sub>CONHNHR]<sub>2</sub> I (R = CHMe<sub>2</sub>, Bu, CH<sub>2</sub>Ph, CH<sub>2</sub>CH<sub>2</sub>Ph, etc.; X = O, S, OCH<sub>2</sub>CH<sub>2</sub>O, RN; n = 1 or 2) were evaluated as monoamine oxidase [9001-66-5] inhibitors. The ED<sub>50</sub> and LD<sub>50</sub> of several compds. in

mice is given. I (R = Bu, X = S, n = 2) [1001-39-4] had the lowest ED50 (2 mg/kg) and TAPHA [I (R = CHMe2, X = S, n = 2)] [1689-03-8] was the least toxic. Structure-activity relations are discussed.

CC 1-3 (Pharmacodynamics)

Section cross-reference(s): 7

IT 999-25-7 999-27-9 999-45-1 1001-28-1 1001-35-0 1001-37-2  
 1001-38-3 1001-39-4 1001-41-8 1001-85-0 1027-27-6 1030-58-6  
 1033-66-5 1050-21-1 1054-53-1 1054-54-2 1057-51-8 1070-55-9  
 1070-56-0 1071-54-1 1071-61-0 1071-68-7 1100-68-1 1100-69-2  
 1102-36-9 1103-69-1 1103-72-6 1105-57-3 1190-60-9 1241-05-0  
 1248-93-7 1250-79-9 1689-03-8 6292-68-8 6292-69-9 75487-18-2  
 75487-19-3 75487-20-6 75487-21-7 75487-22-8 75487-23-9  
 75487-24-0 75487-25-1 75487-26-2 75487-27-3 75487-28-4  
 75487-29-5 75487-30-8 75487-31-9 75487-32-0 75499-16-0


RL: BAC (Biological activity or effector, except adverse); BSU (Biological study, unclassified); BIOL (Biological study)  
 (monoamine oxidase-inhibiting activity of)

IT 75487-20-6


RL: BAC (Biological activity or effector, except adverse); BSU (Biological study, unclassified); BIOL (Biological study)  
 (monoamine oxidase-inhibiting activity of)

RN 75487-20-6 HCAPLUS

CN Propanoic acid, 3-[(2-hydrazino-2-oxoethyl)thio]-, hydrazide (9CI) (CA INDEX NAME)



=> d que  
L15



-O E & SPH NH<sub>2</sub>

and biosensor

NODE ATTRIBUTES:  
DEFAULT MLEVEL IS ATOM  
DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:  
RING(S) ARE ISOLATED OR EMBEDDED  
NUMBER OF NODES IS 11

STEREO ATTRIBUTES: NONE  
L17 166 SEA FILE=REGISTRY SSS FUL L15  
L21 12343 SEA FILE=HCAPLUS ABB=ON PLU=ON BIOSENSORS/CT  
L25 847 SEA FILE=HCAPLUS ABB=ON PLU=ON L17  
L28 1 SEA FILE=HCAPLUS ABB=ON PLU=ON L25 AND (L21 OR BIOSENS? OR  
BIO?(5A)?SENSOR? OR BIOCHIP?)

=> d ibib abs hitind hitstr

L28 ANSWER 1 OF 1 HCAPLUS COPYRIGHT 2004 ACS on STN  
ACCESSION NUMBER: 1997:424393 HCAPLUS  
DOCUMENT NUMBER: 127:158681  
TITLE: Synthesis of a surface-active polyamic acid with  
pendant biological linker molecule for specific  
immobilization of antibodies  
AUTHOR(S): Watson, Hazel; Peltonen, Jouko  
CORPORATE SOURCE: Department of Physical Chemistry, Abo Akademi  
University, Porthansgatan 3-5, FIN-20500, Turku,  
Finland  
SOURCE: → Sensors and Actuators, B: Chemical (1997), B39(1-3),  
261-265  
CODEN: SABCEB; ISSN: 0925-4005  
PUBLISHER: Elsevier  
DOCUMENT TYPE: Journal  
LANGUAGE: English


AB The characterization of a novel amphiphilic polyamic acid derivative intended  
for use in biosensor applications is reported. The  
surface-active polymer, synthesized from a modified diamine and  
1,2,4,5-benzenetetracarboxylic dianhydride (pyromellitic dianhydride),  
comprises a polymer backbone with pendant functional groups, capable of  
specific immobilization of antibodies. The polymer is rendered  
amphiphilic by reaction with octadecylamine in a stoichiometric ratio of  
1:2, i.e., equimolar ratios of acid and amine functionality. The Langmuir  
monolayer is expected to be capable of specifically and effectively  
immobilizing antibody fragments introduced into the subphase. Anal. of  
the chemical structure of the polymer mol. of various mol. wts. and  
characterization of the monolayer are presented.  
CC 9-10 (Biochemical Methods)

Section cross-reference(s): 15  
ST polymer immobilization antibody **biosensor**  
IT **Biosensors**  
Immobilization  
Langmuir films  
(synthesis of surface-active polyamic acid with pendant biol. linker  
mol. for specific immobilization of antibodies)  
IT **193539-15-0P**  
RL: RCT (Reactant); SPN (Synthetic preparation); PREP (Preparation); RACT  
(Reactant or reagent)  
(preparation and reaction with octadecylamine)  
IT **193539-18-3P**  
RL: BSU (Biological study, unclassified); SPN (Synthetic preparation);  
BIOL (Biological study); PREP (Preparation)  
(preparation and use in **biosensor** preparation for antibody  
immobilization)  
IT **193539-13-8P**  
RL: BUU (Biological use, unclassified); RCT (Reactant); SPN (Synthetic  
preparation); BIOL (Biological study); PREP (Preparation); RACT (Reactant  
or reagent); USES (Uses)  
(preparation as linker mol. in **biosensor** preparation for antibody  
immobilization)  
IT **96727-40-1**  
RL: RCT (Reactant); RACT (Reactant or reagent)  
(reaction with pyromellitic anhydride)  
IT **193539-15-0P**  
RL: RCT (Reactant); SPN (Synthetic preparation); PREP (Preparation); RACT  
(Reactant or reagent)  
(preparation and reaction with octadecylamine)  
RN 193539-15-0 HCAPLUS  
CN 1H,3H-Benzo[1,2-c:4,5-c']difuran-1,3,5,7-tetrone, polymer with  
2-[(2,4-diaminophenyl)thio]ethanol (9CI) (CA INDEX NAME)

CM 1

CRN 96727-40-1


CMF C8 H12 N2 O S



CM 2

CRN 89-32-7

CMF C10 H2 O6



IT 193539-18-3P

RL: BSU (Biological study, unclassified); SPN (Synthetic preparation);  
 BIOL (Biological study); PREP (Preparation)  
 (preparation and use in **biosensor** preparation for antibody  
 immobilization)


RN 193539-18-3 HCAPLUS

CN 1H,3H-Benzo[1,2-c:4,5-c']difuran-1,3,5,7-tetrone, polymer with  
 2-[(2,4-diaminophenyl)thio]ethyl 2,5-dihydro-2,5-dioxo-1H-pyrrole-1-  
 hexanoate (9CI) (CA INDEX NAME)

CM 1

CRN 193539-17-2


CMF C18 H23 N3 O4 S



CM 2

CRN 89-32-7

CMF C10 H2 O6




IT 193539-13-8P

RL: BUU (Biological use, unclassified); RCT (Reactant); SPN (Synthetic preparation); BIOL (Biological study); PREP (Preparation); RACT (Reactant or reagent); USES (Uses)

(preparation as linker mol. in **biosensor** preparation for antibody immobilization)

RN 193539-13-8 HCAPLUS

CN 1H-Pyrrole-1-hexanoic acid, 2,5-dihydro-2,5-dioxo-, 2-[(2,4-diaminophenyl)thio]ethyl ester, dihydrochloride (9CI) (CA INDEX NAME)




IT 96727-40-1

RL: RCT (Reactant); RACT (Reactant or reagent)  
(reaction with pyromellitic anhydride)

RN 96727-40-1 HCAPLUS

CN Ethanol, 2-[(2,4-diaminophenyl)thio]- (9CI) (CA INDEX NAME)



REFERENCE COUNT:


15

THERE ARE 15 CITED REFERENCES AVAILABLE FOR THIS



RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L5 STR




NODE ATTRIBUTES:  
 DEFAULT MLEVEL IS ATOM  
 DEFAULT ECLEVEL IS LIMITED

AND Biosensor

GRAPH ATTRIBUTES:  
 RING(S) ARE ISOLATED OR EMBEDDED  
 NUMBER OF NODES IS 6

STEREO ATTRIBUTES: NONE  
 L8 9777 SEA FILE=REGISTRY SSS FUL L5  
 L9 STR



NODE ATTRIBUTES:  
 DEFAULT MLEVEL IS ATOM  
 DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:  
 RING(S) ARE ISOLATED OR EMBEDDED  
 NUMBER OF NODES IS 10

STEREO ATTRIBUTES: NONE  
 L10 17 SEA FILE=REGISTRY SUB=L8 SSS FUL L9  
 L18 10 SEA FILE=HCAPLUS ABB=ON PLU=ON L10  
 L19 2 SEA FILE=HCAPLUS ABB=ON PLU=ON L18 AND (IMMOBIL? OR COAT? OR ATTACH? OR BIOSENS? OR BIOCHIP? OR BIO?(2A)?SENS?)  
 L20 10 SEA FILE=HCAPLUS ABB=ON PLU=ON L18 OR L19  
 L21 12343 SEA FILE=HCAPLUS ABB=ON PLU=ON BIOSENSORS/CT  
 L22 34 SEA FILE=HCAPLUS ABB=ON PLU=ON L21 AND L8  
 L23 7 SEA FILE=HCAPLUS ABB=ON PLU=ON L22 AND ?CARBOHYDR?  
 L24 6 SEA FILE=HCAPLUS ABB=ON PLU=ON L23 NOT L20

=&gt; d 124 ibib abs hitind hitstr 1-6 )

L24 ANSWER 1 OF 6 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 2002:555756 HCAPLUS

DOCUMENT NUMBER: 137:121864

TITLE: Biosensor with covalently attached membrane-spanning proteins

INVENTOR(S): Lakey, Jeremy Hugh

PATENT ASSIGNEE(S): Newcastle University Ventures Limited, UK

SOURCE: PCT Int. Appl., 52 pp.

CODEN: PIXXD2

DOCUMENT TYPE:

Patent

LANGUAGE: English  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:


PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2002057780	A1	20020725	WO 2002-GB222	20020118
W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZM, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM				
RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR, BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG				
EP 1352245	A1	20031015	EP 2002-732154	20020118
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO, MK, CY, AL, TR				
US 2004096895	A1	20040520	US 2003-250682	20031017
PRIORITY APPLN. INFO.:				
			GB 2001-1279	A 20010118
			GB 2001-8947	A 20010410
			WO 2002-GB222	W 20020118
AB	The invention concerns a product comprising: a membrane-spanning protein; a lipid membrane formed from amphiphilic mols. and membrane-spanning protein mols.; and a substrate: characterized in that the membrane protein is directly coupled to the substrate. The invention also provides a method for producing such a product which (i) comprises treating a substrate with a hydrophilic coating agent; (ii) providing at least one membrane-spanning protein; (iii) bringing the protein into contact with the treated substrate under conditions for the coupling of the protein directly to the treated substrate; (iv) adding amphiphilic mols. to the protein-coupled substrate to form a lipid membrane. The product is useful for biosensors, protein arrays and the like.			
IC	ICM G01N033-543			
CC	ICS C12Q001-00; G01N027-333			
IT	9-1 (Biochemical Methods). Section cross-reference(s): 6 Actinobacillus pleuropneumoniae Aeromonas Aeromonas salmonicida Amphiphiles Antibiotics Aquifex aeolicus Bartonella bacilliformis <b>Biosensors</b> Bordetella avium Bordetella pertussis Borrelia burgdorferi Brucella Brucella melitensis Burkholderia cepacia Calymmatobacterium granulomatis Chlamydia Chlamydia trachomatis Chlamydophila pneumoniae Citrobacter freundii Coating materials			

Comamonas acidovorans  
Drugs  
Ectothiorhodospira vacuolata  
Eikenella corrodens  
Enterobacter aerogenes  
Enterobacter cloacae  
Escherichia  
Escherichia coli  
Escherichia fergusonii  
Escherichia hermannii  
Escherichia vulneris  
Genetic methods  
Haemophilus  
Haemophilus actinomycetemcomitans  
Haemophilus ducreyi  
Haemophilus influenzae  
Haemophilus parainfluenzae  
Helicobacter pylori  
Histophilus somni  
Hydrophilicity  
Immobilization, molecular or cellular  
Ionophores  
Klebsiella oxytoca  
Klebsiella pneumoniae  
Legionella pneumophila  
Leptothrix discophora  
Mannheimia haemolytica  
Methylococcus capsulatus  
Moraxella catarrhalis  
Mutation  
Neisseria flavescens  
Neisseria gonorrhoeae  
Neisseria lactamica  
Neisseria meningitidis  
Neisseria polysaccharea  
Neisseria sicca  
Pasteurella multocida  
Pesticides  
Photobacterium profundum  
Pseudomonas aeruginosa  
Pseudomonas fluorescens  
Pseudomonas putida  
Pseudomonas syringae  
Rahnella aquatilis  
Rhodobacter blasticus  
Rhodobacter capsulatus  
Rickettsia prowazeki  
Salmonella typhi  
Salmonella typhimurium  
Self-assembled monolayers  
Serratia marcescens  
Serratia odorifera  
Shigella  
Siler  
Sinorhizobium meliloti  
Siphoviridae  
Sulphydryl group  
Thermotoga maritima  
Treponema pallidum


Treponema phagedenis  
 Vibrio alginolyticus  
 Vibrio cholerae  
 Vibrio parahaemolyticus  
 Xenorhabdus nematophila  
 Yersinia enterocolitica  
 Yersinia pestis  
 $\alpha$ -Helix  
 $\beta$ -Barrel  
 (biosensor with covalently attached membrane-spanning proteins)  
 IT Amino acids, analysis  
 Antibodies  
**Carbohydrates**, analysis  
 DNA  
 Fatty acids, analysis  
 Hormones, animal, analysis  
 Ligands  
 Nucleic acids  
 Peptide nucleic acids  
 RNA  
 Receptors  
 Steroids, analysis  
 cDNA  
 RL: ANT (Analyte); ANST (Analytical study)  
 (biosensor with covalently attached membrane-spanning proteins)  
 IT 52-90-4, Cysteine, properties **87707-01-5** 102281-30-1  
**186133-87-9** 202529-31-5 443299-05-6  
 RL: PRP (Properties)  
 (biosensor with covalently attached membrane-spanning proteins)  
 IT **87707-01-5 186133-87-9**  
 RL: PRP (Properties)  
 (biosensor with covalently attached membrane-spanning proteins)  
 RN 87707-01-5 HCAPLUS  
 CN Hexadecanoic acid, (1R)-1-[3-hydroxy-3-oxido-8-oxo-10-(2-pyridinyldithio)-  
 2,4-dioxo-7-aza-3-phosphadec-1-yl]-1,2-ethanediyl ester (9CI) (CA INDEX  
 NAME)

Absolute stereochemistry.

PAGE 1-A



PAGE 1-B




RN 186133-87-9 HCAPLUS


CN 9-Octadecenoic acid (9Z)-, (1R)-1-[3-hydroxy-3-oxido-8-oxo-10-(2-pyridinyldithio)-2,4-dioxo-7-aza-3-phosphadec-1-yl]-1,2-ethanediyl ester  
(9CI) (CA INDEX NAME)

Absolute stereochemistry.  
Double bond geometry as shown.

PAGE 1-A



PAGE 1-B



REFERENCE COUNT: 7 THERE ARE 7 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L24 ANSWER 2 OF 6 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 2001:731334 HCAPLUS

DOCUMENT NUMBER: 135:269619

TITLE: Colorimetric glycopolythiophene biosensors

INVENTOR(S): Charych, Deborah J.; Myung-Gi-Baek, Deborah J.

PATENT ASSIGNEE(S): The Regents of the University of California, USA

SOURCE: U.S. Pat. Appl. Publ., 38 pp., Cont.-in-part of U.S.  
Ser. No. 461,509.

CODEN: USXXCO

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 11  
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 2001026915	A1	20011004	US 2000-734410	20001211
US 6660484	B2	20031209		
US 6001556	A	19991214	US 1996-592724	19960126
US 6183772	B1	20010206	US 1996-609312	19960301
US 6022748	A	20000208	US 1997-920501	19970829
US 6080423	A	20000627	US 1997-944257	19971006
US 6180135	B1	20010130	US 1997-944323	19971006
US 6468759	B1	20021022	US 1998-33557	19980302
US 6306598	B1	20011023	US 1999-337973	19990621
US 6395561	B1	20020528	US 1999-461509	19991214
US 6485987	B1	20021126	US 2000-500295	20000208

PRIORITY APPLN. INFO.:

US 1992-976697	B2	19921113
US 1993-159927	B2	19931130
US 1994-289384	B2	19940811
US 1994-328237	B2	19941024
US 1995-389475	B2	19950213
US 1998-23898	B3	19950213
US 1996-592724	A3	19960126
US 1996-609312	A2	19960301
US 1997-38383P	P	19970214
US 1997-39749P	P	19970303
US 1997-50496P	P	19970623
US 1997-920501	A3	19970829
US 1997-944323	A2	19971006
US 1998-33557	A2	19980302
US 1998-103344	A2	19980623
US 1999-337973	A2	19990621
US 1999-170190P	P	19991210
US 1999-461509	A2	19991214
US 2000-500295	A2	20000208
US 1992-982189	B2	19921125
US 1997-944257	A3	19971006
US 1998-90266P	P	19980622

AB The present invention relates to methods and compns. for the direct detection of analytes using observable spectral changes in biopolymeric systems. In particular, the present invention allows for the direct colorimetric detection of analytes using color changes that occur in glycopolythiophene polymer systems in response to selective binding of analytes.

IC C12Q001-70; G01N033-554; G01N033-569; A61L002-00; B32B027-04

NCL 435005000

CC 9-1 (Biochemical Methods)

IT Agglutinins and Lectins

Antibodies

Antigens

Biopolymers

Carbohydrates, uses

Cardiolipins

Ceramides

Cerebrosides

Enzymes, uses

Gene

Hormones, animal, uses

Ligands

Lysophosphatidylcholines

Nucleic acids

Phosphatidic acids

Phosphatidylcholines, uses

Phosphatidylethanolamines, uses

Phosphatidylglycerols

Phosphatidylinositols

Phosphatidylserines

Polymers, uses

Polyoxyalkylenes, uses

Proteins, general, uses

Receptors

Sialic acids

Sphingomyelins

Steroids, uses

Trisaccharides

Volatile organic compounds

RL: ARG (Analytical reagent use); DEV (Device component use); ANST  
(Analytical study); USES (Uses)

(colorimetric glycopolythiophene biosensors)

# IT Biosensors

(colorimetric glycopolythiophene; colorimetric glycopolythiophene  
biosensors)

IT 56-12-2, 4-Aminobutanoic acid, reactions 60-32-2, 6-Aminocaproic acid  
107-15-3, Ethylenediamine, reactions 373-44-4, 1,8-Octanediamine  
4781-83-3, 2-Iminothiolane hydrochloride 6066-82-6, n-Hydroxysuccinimide  
6964-21-2, Thiophene-3-acetic acid 7087-68-5, Diisopropylethylamine  
7719-09-7, thionyl chloride 34213-86-0 39001-23-5 58414-52-1  
69492-74-6, Thiophene acetic acid 88829-82-7 125700-67-6, Tbtu  
187146-99-2 187147-00-8 363620-44-4 **363620-46-6**  
363620-47-7

RL: RCT (Reactant); RACT (Reactant or reagent)  
(colorimetric glycopolythiophene biosensors)

IT 81253-66-9P 114815-74-6P 321850-00-4P 321850-01-5P 321850-02-6P  
321850-03-7P 321850-40-2P 321850-42-4P 321850-44-6P 363620-33-1P  
363620-35-3P 363620-37-5P **363620-39-7P**

RL: SPN (Synthetic preparation); PREP (Preparation)  
(colorimetric glycopolythiophene biosensors)

# IT 363620-46-6


RL: RCT (Reactant); RACT (Reactant or reagent)  
(colorimetric glycopolythiophene biosensors)

RN 363620-46-6 HCAPLUS

CN Propanamide, 3-[(2-aminoethyl)thio]-N-[4-[(4,7,8,9-tetra-O-acetyl-3,5-  
dideoxy- $\alpha$ -D-glucopyranosyl)oxy]phenyl]- (9CI) (CA INDEX  
NAME)

Absolute stereochemistry.





IT 363620-39-7P


RL: SPN (Synthetic preparation); PREP (Preparation)  
(colorimetric glycopolythiophene biosensors)

RN 363620-39-7 HCAPLUS


CN 3-Thiopheneacetamide, N-[6-oxo-6-[[2-[[3-oxo-3-[[4-[(4,7,8,9-tetra-O-acetyl-3,5-dideoxy- $\alpha$ -D-glucopyranosyl)oxy]phenyl]amino]propyl]thio]ethyl]amino]hexyl]- (9CI) (CA INDEX NAME)

Absolute stereochemistry. Rotation (+).

PAGE 1-A



PAGE 1-B



ACCESSION NUMBER: 2000:421409 HCAPLUS  
 DOCUMENT NUMBER: 133:40210  
 TITLE: Patterned deposition of antibody-binding proteins for optical diffraction-based biosensors  
 INVENTOR(S): McGrath, Kevin; Kaylor, Rosann M.; Everhart, Dennis S.  
 PATENT ASSIGNEE(S): Kimberly-Clark Worldwide, Inc., USA  
 SOURCE: PCT Int. Appl., 35 pp.  
 CODEN: PIXXD2  
 DOCUMENT TYPE: Patent  
 LANGUAGE: English  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2000036416	A1	20000622	WO 1999-US27727	19991122
W: AE, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CR, CU, CZ, DE, DK, DM, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM RW: GH, GM, KE, LS, MW, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG				
US 2001055754	A1	20011227	US 1998-213713	19981217
US 6579673	B2	20030617		
EP 1141709	A1	20011010	EP 1999-960563	19991122
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO				
AU 762900	B2	20030710	AU 2000-17431	19991122
PRIORITY APPLN. INFO.:			US 1998-213713	A 19981217
			WO 1999-US27727	W 19991122

AB The present invention provides an inexpensive and sensitive device and method for detecting and quantifying analytes present in a medium. The device comprises a metalized film upon which is printed a specific, predetd. pattern of an antibody-binding protein. Upon attachment of a target analyte to select areas of the plastic film upon which the protein is printed, diffraction of transmitted and/or reflected light occurs via the phys. dimensions and defined, precise placement of the analyte. A diffraction image is produced which can be easily seen with the eye or, optionally, with a sensing device. An immunosensor for LH had immobilized protein A printed on a gold/Mylar film. The sensor was reacted with monoclonal antibody to LH  $\beta$ .

IC ICM G01N033-543

ICS G01N021-47

CC 9-1 (Biochemical Methods)

Section cross-reference(s): 2, 15

IT **Biosensors**

(immunol., optical, for LH; patterned deposition of antibody-binding proteins for optical diffraction-based biosensors)

IT **Biosensors**

(immunosensors, optical, for LH; patterned deposition of antibody-binding proteins for optical diffraction-based biosensors)

IT Bacteria (Eubacteria)

**Biosensors**

Candida

Cellophane

Drugs

Drugs of abuse  
Environmental analysis  
Escherichia coli  
Films  
Fungi  
Microspheres  
Optical diffraction  
Salmonella  
Streptococcus pneumoniae  
Virus  
Yeast

(patterned deposition of antibody-binding proteins for optical diffraction-based biosensors)

IT Allergens  
Antibodies

Carbohydrates, analysis  
Carcinoembryonic antigen  
Enzymes, analysis  
Haptens  
Hormones, animal, analysis  
Lipids, analysis  
Nucleic acids  
Polysaccharides, analysis  
Proteins, general, analysis  
Rheumatoid factors

RL: ANT (Analyte); ANST (Analytical study)

(patterned deposition of antibody-binding proteins for optical diffraction-based biosensors)

IT 150244-18-1

RL: RCT (Reactant); RACT (Reactant or reagent)

(thiolation with; patterned deposition of antibody-binding proteins for optical diffraction-based biosensors)


IT 150244-18-1

RL: RCT (Reactant); RACT (Reactant or reagent)

(thiolation with; patterned deposition of antibody-binding proteins for optical diffraction-based biosensors)

RN 150244-18-1 HCAPLUS

CN 3-Pyrrolidinesulfonic acid, 2,5-dioxo-1-[[[1-oxo-6-[[1-oxo-3-(2-pyridinyldithio)propyl]amino]hexyl]oxyl]- (9CI) (CA INDEX NAME)



REFERENCE COUNT:

4

THERE ARE 4 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L24 ANSWER 4 OF 6 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 2000:402097 HCAPLUS

DOCUMENT NUMBER: 133:40221

TITLE: Patterned binding of functionalized microspheres for optical diffraction-based biosensors

INVENTOR(S): Everhart, Dennis S.; Kaylor, Rosann M.; McGrath, Kevin

PATENT ASSIGNEE(S): Kimberly-Clark Worldwide, Inc., USA  
 SOURCE: PCT Int. Appl., 38 pp.  
 CODEN: PIXXD2  
 DOCUMENT TYPE: Patent  
 LANGUAGE: English  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2000034781	A2	20000615	WO 1999-US27671	19991122
WO 2000034781	A3	20000817		
W:				
AE, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CR, CU,				
CZ, DE, DK, DM, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL,				
IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA,				
MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI,				
SK, SL, TJ, TM, TR, TT, TZ, UA, UG, UZ, VN, YU, ZA, ZW, AM, AZ,				
BY, KG, KZ, MD, RU, TJ, TM				
RW:				
GH, GM, KE, LS, MW, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE,				
DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF,				
CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG				
US 6221579	B1	20010424	US 1998-210016	19981211
EP 1137942	A2	20011004	EP 1999-961755	19991122
R:				
AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT,				
IE, SI, LT, LV, FI, RO				
AU 759582	B2	20030417	AU 2000-18271	19991122
US 2001004526	A1	20010621	US 2000-733204	20001208
US 6573040	B2	20030603		

PRIORITY APPLN. INFO.: US 1998-210016 A 19981211  
 WO 1999-US27671 W 19991122

AB The present invention provides an inexpensive and sensitive system and method for detecting analytes present in a medium. The system comprises a diffraction enhancing element, such as functionalized microspheres, which are modified such that they are capable of binding with a target analyte. Addnl., the system comprises a polymer film, which may include a metal coating, upon which is printed a specific, predetd. pattern of analyte-specific receptors. Upon attachment of a target analyte to select areas of the polymer film, either directly or with the diffraction enhancing element, diffraction of transmitted and/or reflected light occurs via the phys. dimensions and defined, precise placement of the analyte. A diffraction image is produced which can be easily seen with the eye or, optionally, with a sensing device. Blue polystyrene particles were conjugated with monoclonal antibody. A gold/Mylar film was blocked with  $\beta$ -casein and then antibody was immobilized in a pattern on the surface. LH sample was mixed with the microparticles and then applied to the sensor. A nitrocellulose disk with a small hole in the center was used to wick away excess fluid and unbound microparticles. A point light source was transmitted through the hole and sensor to create a diffraction image on the other side.

IC ICM G01N033-53

CC 9-1 (Biochemical Methods)

Section cross-reference(s): 2, 3, 15

IT **Biosensors**

(immunol., optical, for LH; patterned binding of functionalized microspheres for optical diffraction-based biosensors)

IT **Biosensors**

(immunosensors, optical, for LH; patterned binding of functionalized microspheres for optical diffraction-based biosensors)

IT Bacteria (Eubacteria)

**Biosensors**

Candida albicans

Cellophane

Chelating agents

Drugs

Drugs of abuse

Environmental analysis

Escherichia coli

Films

Fungi

Microspheres

Optical diffraction

Scanning electron microscopy

Streptococcus pneumoniae

Surfactants

Virus

Yeast

(patterned binding of functionalized microspheres for optical diffraction-based biosensors)

IT Antibodies

Antigens

**Carbohydrates**, analysis

Enzymes, analysis

Hormones, animal, analysis

Lipids, analysis

Nucleic acids

Polysaccharides, analysis

Proteins, general, analysis

RL: ANT (Analyte); ARG (Analytical reagent use); DEV (Device component use); ANST (Analytical study); USES (Uses)

(patterned binding of functionalized microspheres for optical diffraction-based biosensors)

IT 150244-18-1

RL: RCT (Reactant); RACT (Reactant or reagent)

(thiolation with; patterned binding of functionalized microspheres for optical diffraction-based biosensors)


IT 150244-18-1

RL: RCT (Reactant); RACT (Reactant or reagent)

(thiolation with; patterned binding of functionalized microspheres for optical diffraction-based biosensors)

RN 150244-18-1 HCAPLUS

CN 3-Pyrrolidinesulfonic acid, 2,5-dioxo-1-[[1-oxo-6-[[1-oxo-3-(2-pyridinyldithio)propyl]amino]hexyl]oxy]- (9CI) (CA INDEX NAME)



L24 ANSWER 5 OF 6 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1999:405159 HCAPLUS

DOCUMENT NUMBER: 131:41789

TITLE: Optical diffraction biosensor

INVENTOR(S): Everhart, Dennis S.; Jones, Mark L.; Kaylor, Rosann Marie  
 PATENT ASSIGNEE(S): Kimberly-Clark Worldwide, Inc., USA  
 SOURCE: PCT Int. Appl., 39 pp.  
 CODEN: PIXXD2  
 DOCUMENT TYPE: Patent  
 LANGUAGE: English  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9931486	A1	19990624	WO 1998-US26759	19981216
W: AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, UZ, VN, YU, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM				
RW: GH, GM, KE, LS, MW, SD, SZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG				
US 6060256	A	20000509	US 1997-991644	19971216
CA 2309595	AA	19990624	CA 1998-2309595	19981216
AU 9919205	A1	19990705	AU 1999-19205	19981216
AU 760500	B2	20030515		
EP 1040338	A1	20001004	EP 1998-963991	19981216
R: BE, DE, ES, FR, GB, IT, NL, SE				
US 6436651	B1	20020820	US 2000-503554	20000211
PRIORITY APPLN. INFO.:				
			US 1997-991644	A 19971216
			WO 1998-US26759	W 19981216

AB The present invention provides an inexpensive and sensitive device and method for detecting and quantifying analytes present in a medium. The device comprises a metalized film (20) upon which is printed a specific, predetd. pattern of analyte-specific receptors (25). Upon attachment of a target analyte to select areas of the plastic film upon which the receptor is printed, diffraction of transmitted and/or reflected light occurs via the phys. dimensions and defined, precise placement of the analyte. A diffraction image is produced which can be easily seen with the eye or, optionally, with a sensing device.

IC ICM G01N021-47

ICS B41M003-00

CC 9-1 (Biochemical Methods)

IT Bacteria (Eubacteria)

#### Biosensors

Candida

Cellophane

Diapers

Drugs

Escherichia coli

Fungi

Haemophilus influenzae

Hepatitis

Human immunodeficiency virus 1

Human immunodeficiency virus 2

Latex

Neisseria meningitidis

Neoplasm

Optical diffraction

Rous sarcoma virus


Salmonella  
 Streptococcus group A  
 Streptococcus group B  
 Streptococcus pneumoniae  
 Virus  
 Yeast  
 (optical diffraction biosensor)

IT Antibodies  
 Carbohydrates, analysis  
 Carcinoembryonic antigen  
 Enzymes, analysis  
 Glass, analysis  
 Haptens  
 Hormones, animal, analysis  
 Lipids, analysis  
 Nucleic acids  
 Polycarbonates, analysis  
 Polysaccharides, analysis  
 Proteins, general, analysis  
 Rheumatoid factors  
 RL: ANT (Analyte); ANST (Analytical study)  
 (optical diffraction biosensor)

IT 169751-10-4  
 RL: RCT (Reactant); RACT (Reactant or reagent)  
 (optical diffraction biosensor)

IT 169751-10-4  
 RL: RCT (Reactant); RACT (Reactant or reagent)  
 (optical diffraction biosensor)

RN 169751-10-4 HCAPLUS  
 CN 3-Pyrrolidinesulfonic acid, 2,5-dioxo-1-[[1-oxo-6-[[1-oxo-3-(2-pyridinyldithio)propyl]amino]hexyl]oxy]-, monosodium salt (9CI) (CA INDEX NAME)



● Na

REFERENCE COUNT: 6 THERE ARE 6 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L24 ANSWER 6 OF 6 HCAPLUS COPYRIGHT 2004 ACS on STN

ACCESSION NUMBER: 1998:43704 HCAPLUS

DOCUMENT NUMBER: 128:152804

TITLE: Antibody immobilization using heterobifunctional crosslinkers

AUTHOR(S): Shriver-Lake, Lisa C.; Donner, Brian; Edelstein, Rebecca; Breslin, Kristen; Bhatia, Suresh K.; Ligler, France S.

CORPORATE SOURCE: Center for Bio/Molecular Science and Engineering,

Naval Research Laboratory, Washington, DC, 20375-5348,  
USA

SOURCE: Biosensors & Bioelectronics (1997), 12(11), 1101-1106  
CODEN: BBIOE4; ISSN: 0956-5663

PUBLISHER: Elsevier Science Ltd.

DOCUMENT TYPE: Journal

LANGUAGE: English

AB Covalent attachment of functional proteins to a solid support is important for biosensors. One method employs thiol-terminal silanes and hetero-bifunctional crosslinkers such as N-succinimidyl 4-maleimidobutylate (GMBS) to immobilize proteins through amino groups onto glass, silica, silicon or platinum surfaces. In this report, several heterobifunctional crosslinkers are compared to GMBS for their ability to immobilize active antibodies onto glass cover slips at a high d. Antibodies were immobilized at densities of 74-220 ng/cm<sup>2</sup> with high levels of specific antigen binding. **Carbohydrate**-reactive crosslinkers were also compared to GMBS using a fiber optic biosensor to detect fluorescently-labeled antigen. At the concns. tested, the antibodies immobilized with **carbohydrate**-reactive crosslinkers bound more antigen than GMBS immobilized antibodies as indicated by the fluorescence signal.

CC 15-1 (Immunochemistry)  
Section cross-reference(s): 9

IT Immunoglobulins  
RL: ANT (Analyte); ANST (Analytical study)  
(G; amine-reactive and **carbohydrate**-reactive heterobifunctional crosslinkers in immobilization of antibodies to)

IT Antibodies  
RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study); RACT (Reactant or reagent); USES (Uses)  
(amine-reactive and **carbohydrate**-reactive heterobifunctional crosslinkers in immobilization of)

IT **Biosensors**  
Immobilization, biochemical  
(amine-reactive and **carbohydrate**-reactive heterobifunctional crosslinkers in immobilization of antibodies)

IT Glass, biological studies  
RL: BPR (Biological process); BSU (Biological study, unclassified); BIOL (Biological study); PROC (Process)  
(amine-reactive and **carbohydrate**-reactive heterobifunctional crosslinkers in immobilization of antibodies)

IT Crosslinking agents  
(heterobifunctional; amine-reactive and **carbohydrate**-reactive heterobifunctional crosslinkers in immobilization of antibodies)


IT 55750-63-5 58626-38-3 68181-17-9, SPDP 112241-19-7  
**115616-51-8** 157797-94-9 **158913-22-5** 174422-72-1  
RL: RCT (Reactant); RACT (Reactant or reagent)  
(amine-reactive and **carbohydrate**-reactive heterobifunctional crosslinkers in immobilization of antibodies)

IT **115616-51-8 158913-22-5**  
RL: RCT (Reactant); RACT (Reactant or reagent)  
(amine-reactive and **carbohydrate**-reactive heterobifunctional crosslinkers in immobilization of antibodies)


RN 115616-51-8 HCAPLUS

CN Propanoic acid, 3-(2-pyridinyldithio)-, hydrazide (9CI) (CA INDEX NAME)





RN 158913-22-5 HCAPLUS  
CN Propanamide, N-[6-[(2,5-dioxo-1-pyrrolidinyl)oxy]-6-oxohexyl]-3-(2-pyridinyldithio)- (9CI) (CA INDEX NAME)



REFERENCE COUNT: 10 THERE ARE 10 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT